


PDS No. 5421x0	PRODUCT DATA SHEET				Page 1 of 1
Revision 01	EASYstrainer small				
	Greiner Item-No. 5421x0				
Valid for Item-No.:	542100 (sterile)	542120 (sterile)	542140 (sterile)	542170 (sterile)	

1.	Description / Specification	
1.1	Description	<p>EASYstrainer small are cell strainers consisting of a two-component polypropylene housing, with a filtration tissue in four color-codes mesh sizes</p> <ul style="list-style-type: none"> - fits all standard 15 ml tubes and smaller tubes such as 12x75 mm tubes, 5 ml reaction vials or 1.5 ml reaction vials - upper part can be inverted for recovery of retained cells - handle and mantle surface for improved aseptic handling - vent opening for fast filtration - stackable - sterile
1.2	Dimensions	See Customer Drawing
1.3	Volume	Reservoir can hold up to 2 ml volume
1.4	Material / Resin	<p><u>Housing</u>: PP (Polypropylene), free of heavy metal</p> <p><u>Filter</u>: PET (Polyethylene terephthalate), free of heavy metal</p>
1.5	Colour	<p><u>Housing</u>: 542120 (mesh size 20 µm): red 542140 (mesh size 40 µm): green 542170 (mesh size 70 µm): blue 542100 (mesh size 100 µm): yellow</p> <p><u>Filter</u>: natural</p>
1.6	Sterilization	SAL 10 ⁻³
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	Warehouse conditions: +3°C to +40°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	5 years
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Blister	1
3.2	Pieces / Box	50
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality to download

4.	Other Information
	Depending on application or volume used, pre-wetting the filter with 500 µl of buffer to be used is recommended. Discard effluent.

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 00	Date 22 July 2020	Date 23 July 2020	Date 23 July 2020	
Date 28.05.2020	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Illig	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.